

## User Procedure for Fluoromax2 Fluorimeter

1. Before starting to use the fluorimeter, *obtain a UV/Vis absorption spectrum* of the sample to find the absorption properties and the optimum excitation wavelength. The optical absorption of a solution should not be higher than 0.2 at the excitation wavelength and at higher wavelength. Too high absorption will lead to uneven light absorption and internal filter effects. For strongly emitting compounds an absorbance of 0.15 is optimum.

2. Now it's time to use the fluorimeter. Log into the logbook! *NEVER Quit the program and NEVER shut down computer!* If you do so, the monochromators need to be recalibrated. To change from the fluorimeter software to the main Windows surface use the "ALT" + "TAB" keys.

3. Change the directory where your data files will be automatically saved by clicking on:

Main (F4)

    Manage Files

        Change Directory

            \Data\Name\Subfolder hit enter

NB: When inputting information into data fields, leave mouse off of window and toggle between fields using the arrow keys. Always remember to hit enter after inputting a value. To return to main screen, click outside the information field.

4. Insert sample after taking a UV/Vis spectra. Use a clean cell for your sample solution or you will produce useless data.

5. *Emission Scan*

Main (F4)

    Define Experiment

        Emission Scan

            Data Acquisition Parameters

Set Excitation Wavelength from UV/Vis spectra (e.g. Last Absorption Peak )

Set Scan Start approx. 10 nm higher than the excitation wavelength

Leave Experiment, Acquisition Mode, and Auto Zero fields alone

6. Slit widths

ACCS'Y (F7)

For Emission: start with Emiss = 0.2, Excit = 1. Adjust as needed to increase or decrease intensity. For example, if signal is weak, try Emiss = 0.5, Excit = 2

## 7. Start Experiment

Main (F4)

Run Experiment

Input Data File name and Data Title, include slit widths in the title for your record

Leave Experiment File and Experiment Title alone

7. Be sure not to saturate the detector. The counts (Y-axis) should not be higher than  $4 \times 10^6$  in "s" mode.

For a high quality scan, set integration time to 1 sec, or 0.2 for a quick scan

## 8. Excitation Scan

Main (F4)

Define Experiment

Excitation Scan

Data Acquisition Parameters

Set Emission Wavelength from last peak in Emission Scan

Set scan start at a reasonable wavelength (higher than 200 nm) and stop the scan approx. 10 nm BELOW the emission wavelength.

Set acquisition mode to "s/r" to correct for the emission spectrum of the Xe-lamp.

## 9. Slit widths

ACCS'Y (F7)

For the slit width values, use the reverse of what was used for the emission scan. Typically, Emiss = 1, Exc = 0.5.

10. All spectra are automatically saved, if you give a new file name for each new scan. To generate ASCII files to import the data into a graphics program (e.g. IgorPro or EXCEL):

Main (F4)

Other Activities

Output to ASCII

11. **NEVER Quit the program and NEVER shut down computer!** If you do so, the monochromators need to be recalibrated. To change from the fluorimeter software to the main Windows surface use the "ALT" + "TAB" keys. In the Windows surface you can copy your ASCII files to a ZipDisk (100MB) or Floppy Disk for data transfer.